




## Effect of Operating Parameters on Purge & Trap Efficiency

The accepted technique for the analysis of volatile organic pollutants is purge and trap. This technique makes use of a helium carrier gas which is bubbled up through the sample releasing the organics and collecting them onto an adsorbant trap which is then backflushed to the gas chromatograph. The mechanism governing the release of the organics from the water is dependent on the solubility of organics in water and their relative volatility. Research and method development must first optimize the operating parameters in order to achieve efficient, accurate results. Purge time must be long enough to ensure total release of the VOC's from the water and the trap desorption temperature should be hot enough to desorb the organics but not cause thermal breakdown.

The purge time was varied from 1 minute to 10 minutes with a trap desorption temperature of 180°C for 10 minutes. Figures 1 and 2 show the chromatograms of benzene, toluene, o-xylene and ethyl benzene at 1 and 10 minutes represents 200 PPB of each component; after reducing the purge time to 1 minute the recoveries ranged from 11.5% for benzene to 35.9% for o-xylene. A graphic representation of this data can be seen in Figure 2.

In terms of trap desorption temperature, the temperature



was varied from 180°C to 60°C in 30°C increments. No change was observed in organic concentrations until the trap desorption temperature was reduced to 120°C. When the trap temperature was lowered to 60°C, 60% of benzene was recovered while only 17% of ethyl benzene was recovered.

These data reveal the necessity for careful choice of purge and trap parameters in order to obtain accurate and meaningful results.

#### **EQUIPMENT**

##### **PURGE AND TRAP**

330 Sample Concentrator,

purge temperature: ambient  
Internal trap: Tenax TA  
**GAS CHROMATOGRAPHY**  
Varian 3700 equipped with a flame ionization detector  
Column: 50m x 0.25mm SE-54 capillary  
Program: 40°C for 2 min., then 6°C/min. to 100°C  
Carrier: Helium

For more information on this and related applications, we suggest the following readings:

Washall J. W. and Wampler, T P., "Purge and Trap Analysis of Aqueous Samples with Cryofocusing," *Am. Lab.* July 70-74, (1985).

Wampler, T P., Bowe, W. A., Higgins, J. and Levy, E. J., "Systems Approach to Automatic Cryofocusing in Purge and Trap, Headspace, and Pyrolytic Analysis," *Am. Lab.*, August (1985).

© courtesy of CDS Analytical Inc.

#### **Headquarters**

JSB International  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 251 47 53  
F +31 (0) 40 251 47 58

Zoex Europe  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 257 39 72  
F +31 (0) 40 251 47 58

#### **Sales and Service**

Netherlands  
Apolloweg 2B  
8239 DA Lelystad  
T +31 (0) 320 87 00 18  
F +31 (0) 320 87 00 19

Belgium  
Grensstraat 7  
Box 3 1831 Diegem  
T +32 (0) 2 721 92 11  
F +32 (0) 2 720 76 22

Germany  
Max-Planck-Strasse 4  
D-47475 Kamp-Lintfort  
T +49 (0) 28 42 9280 799  
F +49 (0) 28 42 9732 638

UK & Ireland  
Cedar Court,  
Grove Park Business Est.  
White Waltham, Maidenhead  
Berks, SL6 3LW  
T +44 (0) 16 288 220 48  
F +44 (0) 70 394 006 78

[info@go-jsb.com](mailto:info@go-jsb.com)  
[www.go-jsb.com](http://www.go-jsb.com)

